SCHWERT KUNDENINFORMATION

SWORD CUSTOMER INFORMATION

REF 47.803.02 (Hersteller-Art.-Nr.) | Bestellnummer: 6802-72 **C €** 0297 | Klasse **I m** | UMDNS **13-517**

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Dear valued customer!

By purchasing this ratchet, you have acquired a high quality product whose proper handling and use is described as follows.

MANUFACTURER

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Obere Hauptstr. 16-22

78606 Seitingen-Oberflacht, Germany Telephone: +49 (0) 7464 / 98 88 0 +49 (0) 7464 / 98 88 88 Fax: E-mail: info@zepf-dental.com www.zepf-dental.com Internet:

Attention

Please read the information in this leaflet thoroughly. Improper handling and care, as well as misuse, can lead to premature wear or risks for patients and users.

INTENDED PURPOSE

Implants are designed for enossal insertion in the upper and lower jaw. The appropriate transgingival abutments are fixed onto the implants according to the indication and are provided with a prosthetic superstructure. Special matched instruments and the torque ratchet are used for the insertion and prosthetic treatment of the implants. The implant manufacturer's specifications therefore also have to be considered and adhered to under all circumstances.

Attention

Instruments and especially implants may only be used by persons who are specially trained or instructed for this purpose.

CONTRAINDICATION

Insufficient quality of bone mass for insertion of dental implants. Infections or ulcers at the point of insertion in the bone.

UNDESIRABLE SIDE EFFECTS

Due to a patient's metal intolerance, on contact with the instrument, although made of medical stainless steel, a hypersensitivity (allergic) reaction may be triggered. Given such a reaction, the procedure has to be aborted and the necessary steps taken.

Attention

Instruments may only be used by persons who are specially trained or instructed for this purpose.

APPLICATION AND SAFETY INFORMATION

Prior to initial use and before each further use, all instruments have to be completely cleaned, disinfected and

It is very important to inspect every product prior to use for damage and wear, e.g. cracks, breaks or defects. Never use damaged products.

The torque ratchets are delivered unsterile. The products must first be cleaned and sterilized. NON STERILE Compatibility with other products has to be checked in advance.

RECONDITIONING

The ratchet was developed for sterilization in the autoclave and has proven its worth for this kind of sterilization. It firstly has be cleaned and sterilized. The ratchet is classified by RKI as "Critical B". Please observe our general instructions on the preparation of resterilizable instruments in accordance with DIN EN ISO 17664, downloadable from our website www.schwert.com. Compatibility with other products like implants and their abutments has to be checked in advance.

Based on the product design, the materials used and the intended purpose, no limit on the maximum number of executable preparation cycles can be defined.

Attention

The products must be prepared separately from similar or identical instruments to avoid mixing up individual components.

SCHWERT Ratchet

REF 47.803.02 (manufacturer art. no.) | Order number: **6802-72**

MEASUREMENT ACCURACY

The measurement accuracy of the torque ratchet is ensured by the spring used for calibration. Calibration of the ratchet was performed on the basis of the following values:

Characteristic:

Setting at 10 Ncm / tolerance: -1.1 / +1.2 Ncm Setting at 20 Ncm / tolerance: -1.4 / +1.6 Ncm Setting at 30 Ncm / tolerance: -1.5 / +1.8 Ncm Setting at 40 Ncm / tolerance: -1.7 / +2.0 Ncm

The implant manufacturer's specifications always have to be observed.

PREPARATION AND TRANSPORT

The instruments must be stored and transported to the preparation site in a closed container to avoid damage to the instruments and contamination of the environment.

MACHINE CLEANING

Place the ratchet in a strainer bowl on the rack cart or on the inserts of the MIS cart and start the cleaning process

- 01. 1 min. pre-rinsing with cold water
- 02. Empty
- 03. 3 min. pre-rinsing with cold water
- 04. Empty
- 05. 5 min. washing at 55°C (131°F) with 0.5 % alkaline cleaner or 45°C (113°F) with enzymatic cleaner
- 06. Empty
- 07. 3 min. neutralization with warm tap water (> 40°C / 104°F) and neutralizer
- 08. Empty
- 09. 2 min. intermediate rinsing with warm tap water $(> 40^{\circ}C / 104^{\circ}F)$
- 10. Empty

DISINFECTION / MACHINE

Perform machine thermal disinfection taking the national requirements concerning the A0 value (see ISO 15883) into account.

DRYING

Drying the instruments with the drying cycle of the cleaning / disinfection device. If necessary, manual drying can also be achieved with the aid of a lint-free cloth.

MANUAL CLEANING

Prepare a cleaning bath according to the manufacturer's specification.

- 01. Rinse the product under cold tap water (> 40°C / 104°F) until all visible contamination has been removed. Stubborn dirt has to be removed with a soft
- 02. Completely insert the products in the prepared cleaning bath. Comply with the exposure time in the manufacturer's specifications.
- 03. Manually clean the inserted instruments with a soft brush. All surfaces have to be brushed off several times.
- 04. Rinse the products very thoroughly under flowing tap water to remove the cleaning agent without leaving residues.

DISINFECTION / MANUAL

Prepare a disinfection bath according to the disinfection manufacturer's specifications.

Place the instruments in the disinfection bath and observe the suggested exposure time.

Rinse the products very thoroughly with DI water to remove the cleaning agent without leaving residues.

FUNCTIONAL INSPECTION AND PACKAGING

Visual inspection for cleanliness; assembly and functional test according to the operating instructions. Repeat the reconditioning process as necessary until the instrument is visually clean. Standard-compliant packaging of the instruments for sterilization in accordance with ISO 11607 and EN 868.

STERILIZATION

Sterilization of the products with the fractionated pre-vacuum process (according to DIN EN ISO 13060 / DIN EN ISO 17665) under consideration of the respective national

3 pre-vacuum phases with at least 60 mbar pressure. Heating to a sterilization temperature of at least 132°C (270°F); max. 137°C (279°F). Shortest holding time: 3 min. Drying time: at least 10 min.

STORAGE

Store instruments in a dry, clean and dust-free environment. Large temperature fluctuations have to be avoided. This way you avoid the formation of condensate and the resulting corrosion. A shelf life specification or functional

restriction is not given for surgical or dental instruments delivered unsterile following manufacture if they are stored properly.

SCHWERT KUNDENINFORMATION

SCHWERT Ratchet

REF 47.803.02 (manufacturer art. no.) | Order number: 6802-72

REPAIRS

Repairs may only be carried out by the manufacturer, otherwise every quarantee expires for the function.

VALIDATION OF RECONDITIONING

It is the responsibility of the user to validate their procedure appropriately.

BRIEF DESCRIPTION

The ratchet head is lockable to disable the torque limitation. In the case of non-lockable ratchet heads, the torque can be set with the setting screw. The torque is infinitely adjustable from 10 Ncm to 40 Ncm. With the aid of various inserts, the torque ratchet is used with different inserts the manual instrument to prepare the implant bed, to fix the implants in their final position, to screw in the locking screws, gingiva formers and titanium abutment screw.

APPLICATION

The torque ratchet is used with various inserts depending on the application, either with or without torque limitation. Prior to initial use and after every use on the patient, the torque ratchet and its inserts has to be disassembled, cleaned, dried and sterilized after assembly. See "Disinfection, cleaning and sterilizing instruments" in these Instructions for Use.

Important:

After each use, the setting screw has to be unscrewed to the lowest marking to relieve the spring! We recommend annual calibration of the torque ratchet, which can be performed by A. Schweickhardt GmbH & Co. KG. The user is responsible for performing the annual calibration at the manufacturer's site or an accredited test laboratory with suitable testing equipment.

USE WITHOUT TORQUE LIMITATION

Adjust the setting screw to 10 Ncm with several turns.

Pull the ratchet head out of the ratchet handle, turn 90° and click into place on the ratchet handle. The ratchet head can now no longer be bent.

After attaching the appropriate inserts (accessory components for manual surgical instruments) the torque ratchet can be used without the torque function.

USE WITH TORQUE LIMITATION

Set the ratchet head such that it bends under lateral load. Once the ratchet head is locked, the setting screw has to be set to 10 Ncm with a few turns. Then the ratchet head is pulled out of the ratchet handle, turned 90° and clicked into place on the ratchet handle.

Tighten the setting screw until the required torque is visible on the measurement line and coincides with the marking on the ratchet handle. Note: 6 turns on the setting screw corresponds to a difference of 10 Ncm.

Abutments, accessories, etc. can now be screwed in according to the implant manufacturer's specifications until the ratchet head bends. This is a sign for reaching the set torque. Do not continue screwing in the bent state.

EPEE - SPADA

SCHWERT Ratchet

REF 47.803.02 (manufacturer art. no.) | Order number: 6802-72

DISASSEMBLY

Remove inserts from the ratchet head. Completely screw out the setting screw from the ratchet head.

Pull off the ratchet head from the ratchet handle and remove the torque spring.

ASSEMBLY

Push the torque spring into the ratchet handle.

Oil the contact and friction surfaces with medical white oil art. no. 3297.

Push the ratchet head with the threaded rod into the ratchet handle.

Oil the joints and threads with medical white oil art. no. 3297. Wipe off excess oil with a lint-free cloth. Screw in the setting screw onto the ratchet handle down to the lowest marking.

HANDLING

Always treat all instruments with the greatest care while transporting, cleaning, care, sterilization and storage. This particularly applies for sensitive areas.

WARRANTY

Our company exclusively delivers inspected and flawless products to its customers. All our products are designed and manufactured such that the highest quality standards are fulfilled. Liability for products that are modified as compared with the original, misused or improperly treated or used, is excluded.

A. SCHWEICKHARDT GMBH & CO. KG ACCEPTS NO LIABILITY IF THIS CUSTOMER INFORMATION HAS BEEN VERIFIABLY VIOLATED.

EVIDENCE OF DECONTAMINATION

On the basis of statutory regulations and above all to protect our employees, we need a signed "Evidence of decontamination" with every product return. Please ensure that this "Evidence of decontamination" is enclosed filled out and signed with every product return (complaint / repair / other reason for return) and package the product such that there is no risk of injury for our personnel in goods receiving during unpacking.

The form is available for download on our website "http://www.schwert.com/media/downloads/ 080300008_e_Proof_of_Decontamination_AS.pdf".

RECOMMENDED CARE PRODUCTS

3297

Paraffinum Perliquidum, 50 ml (Paraffin Oil / Medical White Oil)

Medical Whiteoil for maintenance of instrument hinges and for use as sharpening oil.

Physically harmless acc. to DAB 10.