
Krekeler Messschieber · Krekeler Caliper

Modern Therapy · Modern Instruments

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1732-10

Inserting Implants gets much easier and effective due the caliper acc. to Prof. Krekeler.

Working is safer, faster and more precise by the combination several functions in one instrument.

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For a maximum stability always the widest possible implant diameter should be chosen. By the exactly 1.5 mm measuring rods, the implants maximum diameter is for direct reading.

Due to the latest scientific research in order to avoid bones resorbtion after the implantological treatment, there should be a wall of at least 1.5 mm of bone material to each side.

1. During the measurement of the interdental gap, it is possible to read directly the maximum implant diameter from the **IN** mark. A correct measurement is performed by applying the rods at the cervix dentis of the neighboring teeth. From the **OUT** mark you may read the the overall measure.

2. The sharpened rods also support the operator in placing the first mark in the bone, to ensure a safe distance to the neighbouring implant / tooth. If there is a neighbour implant it is possible to use one rod as spacer for the minimum distance of the next drill.

3. At the distal ending of the Krekeler Caliper you also find a drilling guide and a 5 and 7 mm graduation to measure the patients alveolar ridge from crestal. This supports an easy and direct evaluation if there is enough bone material left.