
UpDate!

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IDR Sofortige Dentoalveolar Restauration Immediate Dentoalveolar Restoration

Modern Therapy · Modern Instruments

IDR

Immediate Dentoalveolar Restoration

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Fig 1 Initial clinical aspect: there is an inflammation at the surrounding tissue of the left central incisor with root fracture.

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Fig 2 A minimally invasive tooth removal with a periotome and mini-elevators to maintain the gingival architecture.

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Fig 3 Periodontal probing evaluation. Total absence of the buccal bone wall was detected.

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Fig 4 Osteotomy with a palatal approach.

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Fig 5 A Replace Select Tapered TiUnite 5x16 mm (Nobel Biocare, Gothenburg, Sweden) was inserted. The final insertion torque was 50 Ncm.

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Fig 6 Confirmation of the extent of the bone defect.

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Fig 7 The extension of the bone defect was evaluated. Measuring the extent of the bone defect in height and width.

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Figs 8-12 A corticomedullary graft was harvested with a straight chisel at the buccal wall of the wisdom tooth region.

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Figs 13 - 14 The final graft dimensions after remodeling according to the defect size and format.

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Figs 15 - 16 The corticomedullary graft is inserted into the alveolar socket in order to obtain primary stability by juxtaposition.

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Fig 17 Bone marrow is inserted and packed between the implant and the corticomedullary bone graft. Observe the thickness of the new buccal cortical bone.

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Figs 20 - 21 After four months, soft tissue maturation at the gingival architecture is seen. The maintenance of the soft tissue volume can be observed.

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Fig 22 Porcelain crown completed and installed. Aesthetic outcome after 2 years of control.

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Fig 23 CBCT sagittal view 3 years later. There is adequate buccal bone thickness at the region of the implant.

IDR

Immediate Dentoalveolar Restoration

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The selection of the right chisel to remove the graft from the maxillary tuberosity is performed according to the shape of the existing defect in the receiving area and the required bone volume. The most suitable chisel shapes are the straight and the gouge ones of 6, 8 or 10 mm width; as a general rule, the chosen chisel should be 2 mm wider than the graft width to be removed.

These chisels have a millimetered active part with the aim to ease the graft measurement at the moment of removal, allowing the obtention of a bone specimen with suitable length and thickness.

Abbildung Illustration	Artikel-Nr. Item-No.	Artikelbezeichnung Description
	1790-90	IDR Kit mit SCHWERT Protector 266 x 180 x 36 mm REF 2603-03 bestehend aus: IDR Kit with SCHWERT Protector 266 x 180 x 36 mm REF 2603-03 consisting of:
	1675-06	6 mm, Lexer-Mini Osteotom / Osteotome, 18 cm
	1675-08	8 mm, Lexer-Mini Osteotom / Osteotome, 18 cm
	1675-10	10 mm, Lexer-Mini Osteotom / Osteotome, 18 cm
	1676-06	6 mm, Lexer-Mini Hohlmeißel / Gouge, 18 cm
	1676-08	8 mm, Lexer-Mini Hohlmeißel / Gouge, 18 cm
	1676-10	10 mm, Lexer-Mini Hohlmeißel / Gouge, 18 cm
	1697-0Ti	Titan Knochenstopfer Graduation 3/6/9/12/15 mm, Ø 1,0 / 2,0 mm Bone Condenser with exchangeable Titanium tips
	2619-10	Medizinbecher, rostfrei Medicine Cup, stainless steel

CURRICULUM VITAE José Carlos Martins da Rosa

- Graduation course in Dentistry, Federal University of Santa Maria Dental School – UFSM – Santa Maria/RS – Brazil, 1988
- Specialization course in Periodontics, São Paulo Association of Dentists – APCD – Bauru/SP – Brazil, 1991
- Specialization course and Master of Science in Prosthesis, CPO SL Mandic – Center of Dental Research São Leopoldo Mandic – Campinas/SP – Brazil, 2005
- Doctorate Student in Implantology, CPO SL Mandic – Center of Dental Research São Leopoldo Mandic – Campinas/SP – Brazil
- Author of the book: Immediate Dentoalveolar Restoration – Immediate-loaded implants in compromised alveolar sockets, by Santos Publishing (São Paulo, Brazil) in Portuguese and Spanish and by Quintessence Publishing in English
- Private Practice, Caxias do Sul – RS – Brazil